



PFI | **VET**

PFI S4 · PFI M4 · PFI X4


Abberation Balanced
Optimal Stability




MEDICENTUR
— VETERINARY —

Abberation Balanced


PFI S4



PFI M4



PFI X4



Technical specification

Type	Single-piece Aspheric Intraocular Lens for implantation into the capsular bag
Optic	Aspheric – Coma Optimized
PCO protection	360° Special Square Edge (patent pending)
Powers available	+42.0 D
Sterilization	Steam

Material

PFI VET	Hydrophilic – Hydrophobic Copolymer 25% water content with UV Absorber
Refractive index	1.46
ABBE number	58

Geometry

Optic design	Biconvex
Overall length	PFI S4: 12.0 mm PFI M4: 13.0 mm PFI X4: 14.0 mm

Optic diameter	6.5 mm
Haptic thickness	0.4 mm
Haptic angulation	0°

Storage

Temperature	+15 – +35°C
Humidity	15% – 50%
Shelf life	5 years (from sterilization)