


AddOn® refractive



Q functional

ONV function*

Artificial lens function

Distant function

Basis function

Features

One piece, foldable intraocular lens (IOL), biocompatible hydrophilic acrylic (25% water content), as additional IOL in combination with a Basis IOL

- | | |
|----------------------|--------------------------|
| Convex/concave optic | ▶ Cell stop/no IOL touch |
| 4 flex-haptics | ▶ Sulcus fit |
| Non torque design | ▶ Rotational stability |
| Square design | ▶ No iris capture |

Optical power (spherical equivalent):



AddOn® refractive	SEQ	Steps
A4SW00	-10.0 D to +10.0 D	0.25 D

*ONV functions: All functions of the human eye which are required to achieve optimal natural vision (ONV) "without glasses" (see rear page)